# PS1000/PS2000

## Kalibrierung Calibration

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### Calibration PS1000 and PS2000:

**Please note: Calibration would not work when calibration unit is not set the right way. When you want to input the calibration weight in kg you should set the calibration unit to kg in menu. If you need decimal place for the input of your calibration weight (e.g. 310.5 kg), please also check the setting of decimal point in calibration mode. Please check calibration unit and decimal point for calibration mode in Setup-menu:**

When scale is in normal weighting mode, press and hold down UNIT and ON/OFF/ZERO key for 3 seconds until the Setup-menu appears.

→ Display zeigt SETUP

- Use TARE to go through the settings. Please move on to setting P9.

After P9. there is shown the current setting of this parameter.

- 0 means no decimal
- 1 means one decimal digit
- 2, 3, 4 means the number of decimal digits
- 5 means no decimal digit and value multiplied by 10

- Press UNIT to change setting
- Press TARE again to save the setting and move on to next point of setting. Next is setting of calibration unit.

After P10. there is shown the current setting of this parameter.

- 0 means kg
- 1 means lb

- Press UNIT to change setting
- Press TARE again to save the setting and move on to next point of setting.
- Press ON/OFF/ZERO to leave Setup-menu

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### calibration:

1. Switch on the scale and wait for normal weighting mode.
2. Press TARE and ON/OFF/ZERO for three seconds.  
→ Indicator shows CAL-?. Scale is ready for calibration.
3. Press TARE to start calibration. (Exit with ON/OFF/ZERO)
4. Indicator shows CAP. – and after this the full capacity that is set for this scale. Indicator shows d. – and after this the value of division in settings.
5. Indicator shows CAL.P0. Please move away any weight of the scale and press TARE to confirm empty scale.
6. Indicator shows CAL.P1. Please put the first calibration weight onto the scale. Use HOLD/PRINT (next digit) and UNITS (raise digit) to input the value of first calibration weight and use TARE to confirm.
7. Indicator shows CAL.P2. Please put the second calibration weight onto the scale. (The second calibration weight can be same as the first calibration weight, when it is ok for accuracy of calibration). Use HOLD/PRINT (next digit) and UNITS (raise digit) to input the value of first calibration weight and use TARE to confirm.
8. Indicator shows CAL.P0 once again. Please move away any weight of the scale and press TARE to confirm empty scale. The indicator unit will restart and return to normal weighting mode.

**Please note: If indicator shows CAL.Er, an error occurred. Calibration will start again after a short period. It shows CAL.P0 again.**